

Nitrile Examination Gloves



[karex.com.my](https://www.karex.com.my)



NITRILE GLOVES

• 3.2 Mil

• 3.5 Mil

• 3.5 Mil Extended Cuff

- Non-sterile
- Online Chlorinated
- Powder Free
- 100% synthetic nitrile latex
- Size: XS, S, M, L, XL
- Colour: Violet, Skyblue, Black, White

Physical Dimensions

DIMENSIONS	ASTM D6319	EN 455
Length (mm)	min.220 (S) min.230 (M,L)	min. 240
Palm Width (mm)		
S	80 ± 10	80 ± 10
M	95 ± 10	95 ± 10
L	110 ± 10	110 ± 10
Single-Wall Thickness (mm)		
Finger	min. 0.05	N/A
Palm	min. 0.05	N/A

Physical Properties

PROPERTIES	ASTM D6319	EN 455
Tensile Strength (MPa)		
Before Aging	min. 14	N/A
After Aging	min. 14	N/A
Elongation At Break (%)		
Before Aging	min. 500	N/A
After Aging	min. 400	N/A
Median Force At Break (N)		
Before Aging	N/A	min. 6
After Aging	N/A	min. 6

SHELF LIFE • 3 years from the date of manufacturing

STORAGE • Avoid direct sunlight.
• The gloves shall maintain their properties when stored in a dry location.

QUALITY STANDARDS • Conforms to ASTM D6319 and EN 455 standards.
• Manufactured under US FDA 21 CFR Part 820 QSR and ISO 12485:2016.

Get in touch with us at enquiry@karex.com.my

Head Office:
Lot 591 & 594, Persiaran Raja Lumu Pandamaran, Industrial Estate
42000 Port Klang, Selangor, Malaysia.

Tel: +603 3165 6688

Fax: +603 3166 2000

karex